


Bed Exit System (BES)




- ❖ Optimize the lean manpower. Reduced need for care staff to be constantly looking up on the residents. Focus on value-added work.
- ❖ Better monitoring and thus safety for Residents. No need to rely human efforts avoiding any lapses in monitoring and Resident back to safety faster.


Step 1. Pressure pad is placed underneath the mattress



Step 2. User lies on the bed, and the console is ready.



Step 3. When the user leaves the bed, the alarm alerts!



*Scalable:

- Accommodates from 1 bed to 9 beds per room, as well as multiple rooms