

MyGuardian Smart Care System (Bed Exit System)

Smart detection with bigger coverage.

Embrace a more convenient, comprehensive, and scalable Bed Exit solution

Smart Pad:

L/W/T: 150/ 79/ 1.4 cm Sensing: 142/66/0.2 cm

Weight: 1 Kg

Material: Flame-resistant, water-repellent,

splash-proof, antibacterial, and mold-

resistant









Features:

- 1. Bed Exit alert with smart weight detection
- 2. Trend help track health condition
- 3. Sleep tracking
- 4. Durable smart pad
- 5. Minimum workflow interruption
- 6. Wireless communicate via WIFI

In-bed States Enter, Inbed, Side-lying,

Bed Exit Alerts

Multi-stage bed exit & bedside voice alert



Sit up

Exit



24/7 Abnormality Watch

Repositioning, off bed reminder





Leaving Overtime

Used in numerous hospitals and nursing homes across Taiwan since 2017



In Bed Slide Lying

Global Link Tech Pte Ltd

13 Lorong 8 Toa Payoh, #04-02 Braddell Tech, Singapore 319261

: sales@glt.sg : https://glt.sg : +65 91121331

CO/GST REG. NO: 201003084W

MyGuardian Smart Care System (Bed Exit System)



Smart pad

The pad is placed beneath a hospital mattress


Transceiver

Sends real-time data via WIFI; lightweight with voice alert capability.

Server


Nerve center of the system performing data analysis, alerting, managing all connected transceivers.








3
















Global Link Tech Pte Ltd

13 Lorong 8 Toa Payoh, #04-02 Braddell Tech, Singapore 319261 : sales@glt.sg
: https://glt.sg

: +65 91121331
CO/GST REG. NO: 201003084W